COLLAPAT® II

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Bone substitute biomaterial
Based on collagen and hydroxyapatite

Osteoconductive and haemostatic
For repairing bone substance losses

Distributor:

Manufactured by **SYMATESE**
biomatériaux

■ Advantages

- **COLLAPAT® II** is composed of two of the main constituents of the bone collagen and hydroxyapatite.
- **COLLAPAT® II** can be impregnated with an antibiotic solution in the event of septic bone infection.
- **COLLAPAT® II** is osteoconductive:
It enables rapid regeneration tissue growth in the collagen network by means of the hydroxyapatite granules.
- **COLLAPAT® II** is moldable:
Due to its spongy structure, the product takes the consistency of a paste in contact with blood or tissue fluids, enabling optimal adaptation to the shape of the bone surface.
- **COLLAPAT® II** is haemostatic:
When applied in haemorrhagic lesion, **COLLAPAT® II** stops the blood flow in a few minutes.
- **COLLAPAT® II** makes it possible to avoid costly autologous or allogenic transplantation operations.

■ Recommandations for use

COLLAPAT® II must be implanted in the areas in contact with live bone to enable strong bone regeneration.

Shelf-life:

COLLAPAT® II can be stored for three years in its original packaging, from the sterilisation date. The expiry date is printed on each unit box.

■ Presentation

COLLAPAT® II is available in three sizes :

REF. PAT35X6 : 3,5 x 6 x 0,6 cm

REF. PAT7X11 : 7 x 11 x 0,6 cm

REF. PAT1X1 : 1 x 1 x 1 cm

COLLAPAT® II is presented in double packaging, impermeable to light, in unit box form, sterilised with Beta rays.

■ Indications

In orthopaedics:

COLLAPAT® II is used to promote the repair of various types of bone lesions:

- After extracting cortico-spongoid bone fragments.
- After tumour resection.
- When changing an implant.
- In surgical spondylodeses.
- In cases of pseudarthrosis.
- In certain fractures treated by osteosynthesis.

COLLAPAT® II is also used to induce bone substance replacement **in maxillo facial surgery and stomatology**.

Etude de Performance / Performance Study:

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
■ Principle:

The osteoconductive properties of **COLLAPAT® II** were demonstrated by means of a study on rabbits. **COLLAPAT® II** was implanted in a corticocancellous defect 4.2 mm in diameter in the femoral site. The progression of bone regrowth was evaluated on histological sections from specimens included in resin magnified by a factor of 20 from T0 to T3 months.

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T0 : Implantation date

The hydroxyapatite granules are dispersed in the collagen matrix.

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T1 : Primary bone

After 1 month, the collagen has been completely resorbed and the implanted sites are largely colonised with primary neoformed bone tissue, with signs of direct bone-hydroxyapatite particle apposition.

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T3 : Mature bone

After 3 months, the bone reconstruction is complete and is characterised by mature remodeled bone tissue.

Lamelle osseuse et densité ostéocytaire accrue.
Bone lamella and increased osteocytic density.

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- **COLLAPAT® II** enables early bone reconstruction, completed within three months.
- **COLLAPAT® II** is perfectly integrated into the bone without undesirable local reactions.

Fabricant / Manufacturer:

SYMATESE biomatériaux - Z.I. Les Troques - 69630 Chaponost - FRANCE
Tél: +33 (0)4 78 56 72 80 - Fax: +33 (0)4 78 56 00 48